

 **Emperor****EmFilm Plus**
Film Dressing with Absorbent Pad

Featuring a latex-free, semi-permeable polyurethane film membrane which is coated with a hypoallergenic acrylic adhesive, EmFilm Plus Dressings are designed to be breathable with good oxygen and moisture vapour exchange.

Waterproof and impermeable to liquids, bacteria and viruses, the risk of infection is reduced, and the dressing also incorporates an integral central pad consisting of absorbent material.



Facilitates wound healing



Provision of a 'moist wound healing environment'



The patient can shower with the dressing in place



Prevention of contamination by dirt or bacteria, thus a reduced risk of wound infection



Rounded edges and consistent hypoallergenic adhesive help prevent dressing damage and early removal, reducing the number of dressing changes needed.

Call: 0800 008 6646

E-mail: woundcare@rociallehealthcare.com


Visit: www.emperorwound.com



Emperor is a product brand of Rociale Healthcare

Save up to 88% vs other brands


EmFilm Plus 10x15cm
 Cosmopore transparent 9x15cm
 Mepore Film + Pad 9x15cm
 Tegaderm + Pad 9x15cm



Prices Based on DT prices August 2024

Size	Product Code	NHSSC Code	PIP Code	Drug Tariff Price
4cm x 5cm	537698ZD	ELW85027	424-3754	£0.13
5cm x 7cm	537690ZD	ELW1079	424-3762	£0.13
10cm x 10cm	537691ZD	ELW1077	424-3721	£0.15
10cm x 15cm	537692ZD	ELW1068	424-3747	£0.20
10cm x 25cm	537693ZD	ELW1066	424-3739	£0.21
10cm x 30cm	537694ZD	ELW1076	425-3713	£0.26

To view our full portfolio of
Wound Care products, visit
www.emperorwound.com



For **FREE** product samples and more information, please get in touch

Call: 0800 008 6646

E-mail: woundcare@rociallehealthcare.com

Visit: www.emperorwound.com



Emperor is a product brand of Rocialle Healthcare